



MADE IN JAPAN




New

# OPHTHALMIC SCALPEL MICRO FEATHER SLIT KNIFE -RIB&R-

MICRO FEATHER OPHTHALMIC SCALPEL Slit Knife <RIB & R> reduces incision resistance and allows smoother incision than ever before by its exclusively designed RIB-shaped knife and arrowhead-shaped curved edge.

### Improved Operability

Length from blade tip to bent part is shortened from 8mm to 6mm.



Sales unit : 5 pcs per box

Blade width	Plastic Handle	Aluminum Handle
1.4 mm	P-0614RC	No.0614RC
1.6 mm	P-0616RC	No.0616RC
1.8 mm	P-0618RC	No.0618RC
2.0 mm	P-0620RC	No.0620RC
2.2 mm	P-0622RC	No.0622RC

Blade width	Plastic Handle	Aluminum Handle
2.4 mm	P-0624RC	No.0624RC
2.75mm	P-06275RC	No.06275RC
2.8 mm	P-0628RC	No.0628RC
3.0 mm	P-0630RC	No.0630RC

■ Sterilized by Gamma radiation    ■ Bevel up



**FEATHER SAFETY RAZOR CO., LTD.**  
OVERSEAS TRADE DIVISION

3-70, OHYODO MINAMI 3-CHOME, KITA-KU, OSAKA 531-0075, JAPAN  
PHONE: +81-6-6458-1638    FAX: +81-6-6458-1611  
URL [www.feather.co.jp/](http://www.feather.co.jp/)    E-mail [overseas@feather.co.jp](mailto:overseas@feather.co.jp)



9001 Cutting Tool  
13485 Medical Products